

anti- FBXL16 antibody

Product Information

Catalog No.:	FNab03035
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Substrate-recognition component of the SCF(SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

Immunogen information

Immunogen:	F-box and leucine-rich repeat protein 16
Synonyms:	C16orf22, FBL16
Observed MW:	52 kDa
UniprotID :	Q8N461

Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03035(FBXL16 Antibody) at dilution of 1:600