

## anti-Fetuin-A antibody

### Product Information

Catalog No.:	FNab03081
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3E0
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	alpha-2-HS-glycoprotein
Synonyms:	None
Observed MW:	58 kDa
UniprotID :	P02765

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF, IP, IHC
Recommended dilution:	WB: 1:2000-1:10000; IP: 1:1000-1:6000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human liver cancer using FNab03081(AHSG antibody) at dilution of 1:200



IP Result of anti-Fetuin-A (IP:FNab03081, 3ug; Detection:FNab03081 1:3000) with HepG2 cells lysate 6000ug.

HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03081(AHSG antibody) at dilution of 1:5000

