

anti- FKBPL antibody

Product Information

Catalog No.:	FNab03153
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4B10
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in response to X-ray. Regulates p21 protein stability by binding to Hsp90 and p21.

Immunogen information

Immunogen:	FK506 binding protein like
Synonyms:	DIR1, NG7
Observed MW:	42 kDa
Uniprot ID :	Q9UIM3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:5000; IHC: 1:50-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab03153(FKBPL Antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03153(FKBPL antibody) at dilution of 1:2000