

anti- FUT4 antibody

Product Information

Catalog No.:	FNab03248
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May catalyze alpha-1,3 glycosidic linkages involved in the expression of Lewis X/SSEA-1 and VIM-2 antigens.

Immunogen information

Immunogen:	fucosyltransferase 4(alpha(1,3) fucosyltransferase, myeloid-specific)
Synonyms:	ELFT, FCT3A
Observed MW:	95-110 kDa
Uniprot ID :	P22083

Application

Reactivity:	Human
Tested Application:	ELISA
Recommended dilution:	None
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com