
GAGE1 antibody

Product Information

Catalog No.:	FNab03298
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


Antigen, recognized on melanoma by autologous cytolytic T-lymphocytes. Completely silent in normal adult tissues, except testis.

Immunogen information


Immunogen:	G antigen 1
Synonyms:	None
Observed MW:	16 kDa
Uniprot ID :	P0DTW1

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab03298(GAGE1 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03298(GAGE1 antibody) at dilution of 1:500