

anti- GAL3ST3 antibody

Product Information

Catalog No.:	FNab03303
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transfers a sulfate to position 3 of non-reducing beta-galactosyl residues in N-glycans and core2-branched O-glycans. Has high activity towards Gal-beta-1,4-GlcNAc, Gal-beta-1,4(Fuc-alpha-1,3)GlcNAc and lower activity towards Gal-beta-1,3(Fuc-alpha-1,4)GlcNAc.

Immunogen information

Immunogen:	galactose-3-O-sulfotransferase 3
Synonyms:	None
Observed MW:	49-50 kDa
UniprotID :	Q96A11

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03303(GAL3ST3 Antibody) at dilution of 1:600