

anti- Gamma adaptin antibody

Product Information

Catalog No.:	FNab03332
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	adaptor-related protein complex 1, gamma 1 subunit
Synonyms:	ADTG, AP 1 complex subunit gamma 1, AP1G1, CLAPG1, Gamma1 adaptin
Observed MW:	91 kDa
Uniprot ID :	O43747

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab03332 (Gamma adaptin Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03332 (Gamma adaptin antibody) at dilution of 1:400