

anti- GATAD2B antibody

Product Information

Catalog No.:	FNab03366
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transcriptional repressor. Enhances MBD2-mediated repression. Efficient repression requires the presence of GATAD2A. Targets MBD3 to discrete loci in the nucleus. May play a role in synapse development.

Immunogen information

Immunogen:	GATA zinc finger domain containing 2B
Synonyms:	KIAA1150
Observed MW:	69 kDa
UniprotID :	Q8WXI9

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



A549 cells were subjected to SDS PAGE followed by western blot with FNab03366(GATAD2B Antibody) at dilution of 1:1000