

## anti- GHITM antibody

### Product Information

Catalog No.:	FNab03450
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Required for the mitochondrial tubular network and cristae organization. Involved in apoptotic release of cytochrome c.

### Immunogen information

Immunogen:	growth hormone inducible transmembrane protein
Synonyms:	DERP2, MICS1, TMBIM5
Observed MW:	25-27 kDa
UniprotID :	Q9H3K2

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human brain using FNab03450(GHITM antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03450(GHITM antibody) at dilution of 1:400