

anti- GHRHR antibody

Product Information

Catalog No.:	FNab03451
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Receptor for GRF, coupled to G proteins which activate adenylyl cyclase. Stimulates somatotroph cell growth, growth hormone gene transcription and growth hormone secretion.

Immunogen information

Immunogen:	growth hormone releasing hormone receptor
Synonyms:	GHRFR, GHRH receptor, GHRHR, GHRHRpsv, GRF receptor, GRFR
Observed MW:	48 kDa
Uniprot ID :	Q02643

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03451(GHRHR antibody) at dilution of 1:500