

anti- GKN1 antibody

Product Information

Catalog No.:	FNab03473
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6F7
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Has mitogenic activity and may be involved in maintaining the integrity of the gastric mucosal epithelium.

Immunogen information

Immunogen:	gastrokine 1
Synonyms:	AMP18, CA11
Observed MW:	22 kDa
UniprotID :	Q9NS71

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:5000-1:50000; IHC: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab03473(GKN1 Antibody) at dilution of 1:200



human stomach tissue were subjected to SDS PAGE followed by western blot with FNab03473(GKN1 Antibody) at dilution of 1:20000