

anti- GOSR2/Membrin antibody

Product Information

Catalog No.: FNab03566

Size: 100μg Form: liquid

Purification: Protein A+G purification

Purity: ≥95% as determined by SDS-PAGE

Host: Mouse

Clonality: monoclonal

Clone ID: 8G6
IsoType: IgG2a

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

This gene encodes a trafficking membrane protein which transports proteins among the medialand trans-Golgi compartments. Due to its chromosomal location and trafficking function, this gene may be involved in familial essential hypertension.

Immunogen information

Immunogen: golgi SNAP receptor complex member 2

Synonyms: 27 kDa Golgi SNARE protein, Bos1, GOSR2, GS27, Membrin

Observed MW: 25 kDa Uniprot ID: 014653

Application

Reactivity: Human

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human liver slide using FNab03566(GOSR2 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03566(GOSR2 antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>