

anti- GPR3 antibody

Product Information

Catalog No.:	FNab03604
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Orphan receptor with constitutive G(s) signaling activity that activate cyclic AMP. Has a potential role in modulating a number of brain functions, including behavioral responses to stress(By similarity), amyloid-beta peptide generation in neurons and neurite outgrowth(By similarity). Maintains also meiotic arrest in oocytes(By similarity).

Immunogen information

Immunogen:	G protein-coupled receptor 3
Synonyms:	ACCA
Observed MW:	None
Uniprot ID :	P46089

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab03604(GPR3 antibody) at dilution of 1:50