

anti- GPR50 antibody

Product Information

Catalog No.:	FNab03607
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Does not bind melatonin.

Immunogen information

Immunogen:	G protein-coupled receptor 50
Synonyms:	G protein coupled receptor 50, GPR50, H9, Mel1c, Melatonin related receptor
Observed MW:	67 kDa
Uniprot ID :	Q13585

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IF, IHC, FC
Recommended dilution:	IHC: 1:20-1:200; IF: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab03607(GPR50Antibody) at dilution of 1:50



Immunofluorescent analysis of SH-SY5Y cells using FNab03607(GPR50Antibody) at dilution of 1:100