

## anti- GPRC5A antibody

### Product Information

Catalog No.:	FNab03614
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	G protein-coupled receptor, family C, group 5, member A
Synonyms:	GPCR5A, GPRC5A, RAI3, RAIG 1, RAIG1
Observed MW:	40-45 kDa
Uniprot ID :	Q8NFI5

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200 - 1:2000; IHC: 1:50 - 1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab03614 (GPRC5A, RAI3 Antibody) at dilution of 1:200



A375 cells were subjected to SDS PAGE followed by western blot with FNab03614 (GPRC5A, RAI3 Antibody) at dilution of 1:2000