

anti-Granulin antibody

Product Information

Catalog No.:	FNab03631
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	granulin
Synonyms:	acroganin, GEP, GP88, granulin, Granulin A Granulin 5, Granulin B Granulin 4, Granulin C Granulin 6, Granulin D Granulin 7, Granulin E, Granulin F Granulin 3, Granulin G Granulin 2, Granulins, GRN, PCDGF, PEPI, PGRN, Proepithelin
Observed MW:	78 kDa
UniprotID :	P28799

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human prostate cancer slide using FNab03631(PCDGF,GRN Antibody) at dilution of 1:50



Hela cells were subjected to SDS PAGE followed by western blot with FNab03631(PCDGF,GRN antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com