

anti- GSNOR,ADH5 antibody

Product Information

Catalog No.:	FNab03679
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3E7
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	alcohol dehydrogenase 5(class III), chi polypeptide
Synonyms:	ADH 3, ADH5, ADHX, Alcohol dehydrogenase 5, Alcohol dehydrogenase class 3, FALDH, FDH, GSH FDH, GSNOR, GSNOR,ADH5
Observed MW:	40 kDa
UniprotID :	P11766

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:1000-1:8000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275


Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab03679(GSNOR,ADH5 Antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells using FNab03679 (GSNOR,ADH5 antibody) at dilution of 1:50 and Rhodamine-Goat anti-Mouse IgG



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab03679(GSNOR,ADH5 Antibody) at dilution of 1:4000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com