

anti- Hexokinase 2 antibody

Product Information

Catalog No.:	FNab03848
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	hexokinase 2
Synonyms:	Hexokinase 2, Hexokinase type II, HK II, HK2, HKII, HXK2, Muscle form hexokinase
Observed MW:	102kd
Uniprot ID :	P52789

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:2000-1:10000; IP: 1:500-1:2000; IHC: 1:100-1:600
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889


www.fn-test.com



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab03848(HK2 antibody) at dilution of 1:50



IP Result of anti-Hexokinase 2 (IP: FNab03848, 3ug; Detection: FNab03848 1:1000) with HeLa cells lysate 2000ug.



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03848(HK2 antibody) at dilution of 1:1000