

anti- HIF1AN antibody

Product Information

Catalog No.:	FNab03860
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Immunogen information

Immunogen:	hypoxia inducible factor 1, alpha subunit inhibitor
Synonyms:	FIH1
Observed MW:	40 kDa
UniprotID :	Q9NWT6

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200; IP : 1:50 - 1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded rat brain using FNab03860(HIF1AN antibody) at dilution of 1:100



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab03860(HIF1AN antibody) at dilution of 1:1000