

1

anti- HMOX1 antibody

Product Information

Catalog No.:	FNab03937
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Heme oxygenase cleaves the heme ring at the alpha methene bridge to form biliverdin. Biliverdin is subsequently converted to bilirubin by biliverdin reductase. Under physiological conditions, the activity of heme oxygenase is highest in the spleen, where senescent erythrocytes are sequestrated and destroyed. Exhibits cytoprotective effects since excess of free heme sensitizes cells to undergo apoptosis.

Immunogen information

Immunogen:	hemeoxygenase(decycling) 1
Synonyms:	Heme oxygenase 1, HMOX1, HO, HO 1, HO1, HSP32
Observed MW:	33 kDa
UniprotID :	P09601

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP, IF
Recommended dilution:WB: 1:500-1:1000; IP: 1:500-1:2000; IHC: 1:100-1:400; IF: 1:20-1:200	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com




Image:



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab03937(HO-1 Antibody) at dilution of 1:200

Immunofluorescent analysis of (-20°C Ethanol) fixed HT-1080 cells using FNab03937 (HO-1 antibody) at dilution of 1:50 and Alexa Fluor 488congugated Goat Anti-Rabbit IgG(H+L)




IP Result of anti-HMOX1 (IP:FNab03937, 3ug; Detection:FNab03937 1:1000) with HeLa cells lysate 3000ug.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027





mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab03937(HMOX1 antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>