

# anti- HMOX2 antibody

#### **Product Information**

Catalog No.: FNab03938

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Heme oxygenase cleaves the heme ring at the alpha methene bridge to form biliverdin. Biliverdin is subsequently converted to bilirubin by biliverdin reductase. Under physiological conditions, the activity of heme oxygenase is highest in the spleen, where senescent erythrocytes are sequestrated and destroyed. Heme oxygenase 2 could be implicated in the production of carbon monoxide in brain where it could act as a neurotransmitter.

#### **Immunogen information**

Immunogen: hemeoxygenase(decycling) 2

Synonyms: hemeoxygenase(decycling) 2, Heme oxygenase 2, HMOX2, HO 2, HO2

Observed MW: 36 kDa UniprotID: P30519

#### **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IHC, IF, WB, IP

Recommended dilution: WB: 1:500-1:1000; IP: 1:500-1:1000 IHC: 1:50-1:500; IF: 1:50-1:500

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com


1




# Image:



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using FNab03938(HMOX2 Antibody) at dilution of 1:200.



Immunofluorescent analysis of Hela cells, using HMOX2 antibody FNab03938 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye)



IP Result of anti-HMOX2 (IP:FNab03938, 4ug; Detection:FNab03938 1:500) with Raji cells lysate 2000ug.


# 2

# Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03938(HMOX2 Antibody) at dilution of 1:600

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com