

anti- ICAM3 antibody

Product Information

Catalog No.:	FNab04106
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

ICAM proteins are ligands for the leukocyte adhesion protein LFA-1(integrin alpha-L/beta-2)(PubMed:1448173). ICAM3 is also a ligand for integrin alpha-D/beta-2. In association with integrin alpha-L/beta-2, contributes to apoptotic neutrophil phagocytosis by macrophages(PubMed:23775590).

Immunogen information

Immunogen:	intercellular adhesion molecule 3
Synonyms:	None
Observed MW:	110-160 kDa, 60 kDa
UniprotID :	P32942

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:WB:	1:500-1:2000; IHC: 1:20-1:200


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human colon using FNab04106(ICAM3 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab04106(ICAM3 antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com