

anti- IκB alpha antibody

Product Information

Catalog No.:	FNab04198
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, alpha
Synonyms:	I kappa B alpha, IkappaBalpha, IκB alpha, IKBA, Iκβ-α, MAD 3, MAD3, NF kappa B inhibitor alpha, NFKBI, NFKBIA
Observed MW:	39 kDa
Uniprot ID :	P25963

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:50-1:100
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab04198(IKBA Antibody) at dilution of 1:100



A549 cells were subjected to SDS PAGE followed by western blot with FNab04198(IKBA antibody) at dilution of 1:1000