

anti- IL12A antibody

Product Information

Catalog No.:	FNab04218
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3C7
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated Killer cells, and stimulate the production of IFN-gamma by resting PBMC.

Immunogen information

Immunogen:	interleukin 12A(natural killer cell stimulatory factor 1, cytotoxic lymphocyte maturation factor 1, p35)
Synonyms:	CLMF, CLMF p35, IL 12 subunit p35, IL 12A, IL12A, Interleukin 12 subunit alpha, NFSK, NKSF1, P35
Observed MW:	25kd
Uniprot ID :	P29459

Application

Reactivity:	Human
Tested Application:	ELISA, WB,IHC
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:50-1:500

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889

www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab04218(IL12A Antibody) at dilution of 1:200



Recombinant protein were subjected to SDS PAGE followed by western blot with FNab04218(IL12A Antibody) at dilution of 1:2000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com