

## anti- IL-16 antibody

### Product Information

Catalog No.:	FNab04223
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	5A9
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	interleukin 16(lymphocyte chemoattractant factor)
Synonyms:	HsT19289, IL 16, IL16, LCF, prIL 16, Pro interleukin 16
Observed MW:	37 kDa
Uniprot ID :	Q14005

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:5000-1:20000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Recombinant protein were subjected to SDS  
PAGE followed by western blot with  
FNab04223(IL-16 antibody) at dilution of 1:10000

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)