

anti- IL20RB antibody

Product Information

Catalog No.:	FNab04249
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

The IL20RA/IL20RB dimer is a receptor for IL19, IL20 and IL24. The IL22RA1/IL20RB dimer is a receptor for IL20 and IL24.

Immunogen information

Immunogen:	interleukin 20 receptor beta
Synonyms:	DIRS1
Observed MW:	45kd,66-70kd
UniprotID :	Q6UXL0

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab04249(IL20RB antibody) at dilution of 1:500