

anti- IL-28A antibody

Product Information

Catalog No.:	FNab04262
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1H10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	interleukin 28A(interferon, lambda 2)
Synonyms:	Cytokine Zcyto20, IFN lambda 2, IFNL2, IL 28A, IL28A, IL-28A, Interferon lambda 2, Interleukin 28A, ZCYTO20
Observed MW:	34 kDa
Uniprot ID :	Q8IZJ0

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:5000-1:50000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab04262(IL28A Antibody) at dilution of 1:50



Recombinant protein were subjected to SDS PAGE followed by western blot with FNab04262(IL28A Antibody) at dilution of 1:5000