

anti- IL-29 antibody

Product Information

Catalog No.:	FNab04265
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6G0
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	interleukin 29(interferon, lambda 1)
Synonyms:	None
Observed MW:	22kDa
UniprotID :	Q8IU54

Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:5000-1:20000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



GST tagged recombinant IL29 protein were subjected to SDS PAGE followed by western blot with FNab04265(IL-29 antibody at dilution of 1:10000