

anti- IMP3 antibody

Product Information

Catalog No.:	FNab04297
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4D1
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Component of the 60-80S U3 small nucleolar ribonucleoprotein(U3 snoRNP). Required for the early cleavages during pre-18S ribosomal RNA processing.

Immunogen information

Immunogen:	IMP3, U3 small nucleolar ribonucleoprotein, homolog(yeast)
Synonyms:	C15orf12, MRPS4
Observed MW:	22 kDa
UniprotID :	Q9NV31

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:50-1:500; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab04297(IMP3 Antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04297(IMP3 Antibody) at dilution of 1:2000