

anti- INF2 antibody

Product Information

Catalog No.:	FNab04307
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Severs actin filaments and accelerates their polymerization and depolymerization.

Immunogen information

Immunogen:	inverted formin, FH2 and WH2 domain containing
Synonyms:	C14orf151, C14orf173
Observed MW:	180-200 kDa
UniprotID :	Q27J81

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney using FNab04307(INF2 antibody) at dilution of 1:100



A431 cells were subjected to SDS PAGE followed by western blot with FNab04307(INF2 antibody) at dilution of 1:500