

anti- Inhibin beta A-Specific antibody

Product Information

Catalog No.:	FNab04319
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	inhibin, beta A
Synonyms:	Activin beta A chain, EDF, FRP, INHBA, INHBA-Specific, Inhibin beta A chain, inhibin, beta A
Observed MW:	47 kDa
Uniprot ID :	P08476

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



IP Result of anti-Inhibin beta A-Specific
(IP:FNab04319, 4ug; Detection:FNab04319 1:500)
with rat brain tissue lysate 4000ug.



rat brain tissue were subjected to SDS PAGE
followed by western blot with FNab04319(Inhibin
beta A-Specific antibody) at dilution of 1:500