

## anti- Integrin alpha-3 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab04345   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Protein A+G purification  |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Mouse   |
| Clonality:    | monoclonal  |
| Clone ID:     | 4F6   |
| IsoType:      | IgG2a   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

ITGA3, also named as MSK18, CD49c, GAPB3 and FRP-2, belongs to the integrin alpha chain family. It can be cleaved into heavy chain(840aa) and light chain(176aa). This antibody can recognize the heavy chain.

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | integrin, alpha 3   |
| Synonyms:    | CD49C, FRP 2, Galactoprotein B3, GAP B3, GAPB3, Integrin alpha 3, ITGA3, MSK18, VCA 2, VL3A, VLA 3 subunit alpha, VLA3a |
| Observed MW: | 120 kDa   |
| UniprotID :  | P26006  |

### Application

|                          |              |
|--------------------------|--------------|
| Reactivity:              | Human, Mouse |
| Tested Application:      | ELISA, WB    |
| Recommended dilution:WB: | 1:500-1:2000 |


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Image:



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab04345(ITGA3 antibody) at dilution of 1:1000

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)