

anti- KATNAL1 antibody

Product Information

Catalog No.:	FNab04470
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Regulates microtubule dynamics in Sertoli cells, a process that is essential for spermiogenesis and male fertility. Severs microtubules in an ATP-dependent manner, promoting rapid reorganization of cellular microtubule arrays (By similarity).

Immunogen information

Immunogen:	katanin p60 subunit A-like 1
Synonyms:	None
Observed MW:	50-55 kDa
UniprotID :	Q9BW62

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-KATNAL1 (IP:FNab04470, 3ug;
Detection:FNab04470 1:500) with HeLa cells
lysate 3000ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab04470(KATNAL1
Antibody) at dilution of 1:300