

## anti-Kindlin 2 antibody

### Product Information

Catalog No.:	FNab04576
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	fermitin family homolog 2(Drosophila)
Synonyms:	None
Observed MW:	78 kDa
UniprotID :	Q96AC1

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500-1:3000; IP: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab04576(FERMT2 antibody) at dilution of 1:50



IP Result of anti-FERMT2 (IP:FNab04576, 3ug; Detection:FNab04576 1:1000) with A549 cells lysate 3300ug.

A549 cells were subjected to SDS PAGE followed by western blot with FNab04576(FERMT2 antibody) at dilution of 1:500

