

anti- KPTN antibody

Product Information

Catalog No.:	FNab04639
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Necessary for normal neuromorphogenesis. May be involved in actin dynamics. May play a role in producing the sensory apparatus in hair cells. May play a role in actin rearrangements that accompany platelet activation and stereocilia formation.

Immunogen information

Immunogen:	kaptin(actin binding protein)
Synonyms:	None
Observed MW:	48 kDa
UniprotID :	Q9Y664

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney using FNab04639(KPTN antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04639(KPTN antibody) at dilution of 1:400