

anti- KRAP, SSFA2 antibody

Product Information

Catalog No.:	FNab04640
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	sperm specific antigen 2
Synonyms:	None
Observed MW:	180 kDa
Uniprot ID :	P28290

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IP: 1:500-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab04640(KRAP, SSFA2 Antibody) at dilution of 1:200



IP Result of anti-KRAP, SSFA2 (IP:FNab04640, 4ug; Detection:FNab04640 1:500) with mouse testis tissue lysate 4000ug.

rat testis tissue were subjected to SDS PAGE followed by western blot with FNab04640(KRAP, SSFA2 antibody) at dilution of 1:600

