

## anti- KRT12 antibody

### Product Information

Catalog No.:	FNab04650
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May play a unique role in maintaining the normal corneal epithelial function. Together with KRT3, essential for the maintenance of corneal epithelium integrity(By similarity).

### Immunogen information

Immunogen:	keratin 12
Synonyms:	CK 12, Cytokeratin 12, K12, keratin 12, KRT12
Observed MW:	50-55 kDa
Uniprot ID :	Q99456

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



mouse eye tissue were subjected to SDS PAGE followed by western blot with FNab04650(KRT12 antibody) at dilution of 1:800