

## anti- LAMP3 antibody

### Product Information

Catalog No.:	FNab04687
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

### Background

May play a role in dendritic cell function and in adaptive immunity.

### Immunogen information

Immunogen:	lysosomal-associated membrane protein 3
Synonyms:	CD208, DC LAMP, DCLAMP, LAMP, LAMP 3, LAMP3, Protein TSC403, TSC403
Observed MW:	70 kDa
Uniprot ID :	Q9UQV4

### Application

Reactivity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human lung using FNab04687(LAMP3 antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab04687(LAMP3 antibody) at dilution of 1:1000