

## anti- LGALS3BP antibody

### Product Information

Catalog No.:	FNab04758
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Promotes integrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells.

### Immunogen information

Immunogen:	lectin, galactoside-binding, soluble, 3 binding protein
Synonyms:	90K, BTBD17B, Galectin 3 binding protein, LGALS3BP, M2BP, Mac 2 binding protein, MAC 2 BP, MAC2BP, Tumor associated antigen 90K
Observed MW:	65-70 kDa
UniprotID :	Q08380

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human breast cancer slide using FNab04758(LGALS3BP Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab04758(LGALS3BP antibody) at dilution of 1:1000