

anti- LPO antibody

Product Information

Catalog No.:	FNab04834
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Antimicrobial agent which utilizes hydrogen peroxide and thiocyanate(SCN) to generate the antimicrobial substance hypothiocyanous acid(HOSCN)(By similarity). May contribute to airway host defense against infection.

Immunogen information

Immunogen:	lactoperoxidase
Synonyms:	SAPX
Observed MW:	80 kDa
UniprotID :	P22079

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab04834(LPO Antibody) at dilution of 1:50



human saliva tissue were subjected to SDS PAGE followed by western blot with FNab04834(LPO Antibody) at dilution of 1:1000