

## anti- MARCH8 antibody

### Product Information

Catalog No.:	FNab05004
Size:	100 $\mu$ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

E3 ubiquitin-protein ligase that mediates ubiquitination of CD86 and MHC class II proteins, such as HLA-DR alpha and beta, and promotes their subsequent endocytosis and sorting to lysosomes via multivesicular bodies. May also promote ubiquitination and endocytosis of TFRC and FAS.

### Immunogen information

Immunogen:	membrane-associated ring finger(C3HC4) 8
Synonyms:	MIR, RNF178
Observed MW:	33kd
Uniprot ID :	Q5T0T0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:200-1:1000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



IP Result of anti-MARCH8 (IP:FNab05004, 3ug;  
Detection:FNab05004 1:400) with mouse lung  
tissue lysate 4000ug.



mouse testis tissue were subjected to SDS PAGE  
followed by western blot with  
FNab05004(MARCH8 antibody) at dilution of  
1:800