

anti- MIP-3-beta antibody

Product Information

Catalog No.:	FNab05194
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	chemokine(C-C motif) ligand 19
Synonyms:	Beta chemokine exodus 3, C C motif chemokine 19, CCL19, CK beta 11, CKb11, EBI1 ligand chemokine, ELC, MIP 3 beta, MIP 3b, MIP3B, MIP-3-beta, SCYA19, Small inducible cytokine A19
Observed MW:	12 kDa
Uniprot ID :	Q99731

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500; IF: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using FNab05194(MIP-3-beta Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



A549 cells were subjected to SDS PAGE followed by western blot with FNab05194(MIP-3-beta antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com