

anti- MKS3 antibody

Product Information

Catalog No.:	FNab05208
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	transmembrane protein 67
Synonyms:	None
Observed MW:	None
UniprotID :	Q5HYA8

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IF, IHC
Recommended dilution:	IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded
human kidney tissue slide using
FNab05208(MKS3 antibody)at dilution of 1:50