

anti- MMP21 antibody

Product Information

Catalog No.:	FNab05240
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays a specialized role in the generation of left-right asymmetry during embryogenesis. May act as a negative regulator of the NOTCH-signaling pathway(PubMed:26429889, PubMed:26437028). Cleaves alpha-1-antitrypsin(PubMed:12617721).

Immunogen information

Immunogen:	matrix metalloproteinase 21
Synonyms:	None
Observed MW:	67 kDa
UniprotID :	Q8N119

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human skin cancer tissue slide using FNab05240(MMP21 Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab05240(MMP21 Antibody) at dilution of 1:600