

anti- MRPL24 antibody

Product Information

Catalog No.:	FNab05323
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	mitochondrial ribosomal protein L24
Synonyms:	None
Observed MW:	25 kDa
UniprotID :	Q96A35

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

IP Result of anti-MRPL24 (IP:FNab05323, 4ug;
Detection:FNab05323 1:1000) with HEK-293 cells
lysate 1480ug.



mouse liver tissue were subjected to SDS PAGE
followed by western blot with
FNab05323(MRPL24 antibody) at dilution of
1:500