

# anti- MYBPC3 antibody

## **Product Information**

Catalog No.:	FNab05455	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	$\geq$ 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

### Background

Thick filament-associated protein located in the crossbridge region of vertebrate striated muscle a bands. In vitro it binds MHC, F-actin and native thin filaments, and modifies the activity of actin-activated myosin ATPase. It may modulate muscle contraction or may play a more structural role.

#### **Immunogen information**

Immunogen:	myosin binding protein C, cardiac
Synonyms:	MYBPC3, DKFZp779E1762, CMH4, FHC, Cardiac MyBP C
Observed MW:	141-150kd
Uniprot ID :	Q14896

## Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB		
Recommended dilution: WB: 1:500-1:2000			
Image:			

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan,<br/>Hubei, China(430206)Tel :( 0086)027-87384275Fax: (0086)027-87800889www.fn-test.com

1





Various tissues were subjected to SDS PAGE followed by western blot with FNab05455(MYBPC3 antibody) at dilution of 1:500

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :( 0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com