

anti- Myoglobin antibody

Product Information

Catalog No.:	FNab05515
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4C8
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	myoglobin
Synonyms:	MB, MYO, myoglobin
Observed MW:	18 kDa
Uniprot ID :	P02144

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:5000-1:50000; IHC: 1:50-1:300
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab05515(MB Antibody) at dilution of 1:200



human heart tissue were subjected to SDS PAGE followed by western blot with FNab05515(MB Antibody) at dilution of 1:40000