

## anti- NAT1 antibody

### Product Information

Catalog No.:	FNab05561
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Participates in the detoxification of a plethora of hydrazine and arylamine drugs. Catalyzes the N-or O-acetylation of various arylamine and heterocyclic amine substrates and is able to bioactivate several known carcinogens.

### Immunogen information

Immunogen:	N-acetyltransferase 1(arylamine N-acetyltransferase)
Synonyms:	AAC1
Observed MW:	32-34 kDa
UniprotID :	P18440

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human breast cancer using FNab05561(NAT1 antibody) at dilution of 1:50



A375 cells were subjected to SDS PAGE followed by western blot with FNab05561(NAT1 antibody) at dilution of 1:300