

## anti- NKG2A antibody

### Product Information

Catalog No.:	FNab05744
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	killer cell lectin-like receptor subfamily C, member 1
Synonyms:	CD159A, KLRC1, NK cell receptor A, NKG2, NKG2A
Observed MW:	33kd
UniprotID :	P26715

### Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:1000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05744(NKG2A antibody) at dilution of 1:800