

anti- NOB1 antibody

Product Information

Catalog No.:	FNab05779
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8G0
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May play a role in mRNA degradation.

Immunogen information

Immunogen:	NIN1/RPN12 binding protein 1 homolog(<i>S. cerevisiae</i>)
Synonyms:	ART4, NOB1P, PSMD8BP1
Observed MW:	47kd
UniprotID :	Q9ULX3

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:1000-1:6000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab05779(NOB1 Antibody) at dilution of 1:200



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05779(NOB1 Antibody) at dilution of 1:3000